


ULTRAVIOLET-VISIBLE

SPECTROSCOPY


WHAT HAPPENS WHEN A MOLECULE ABSORBS UV-VISIBLE RADIATION?

- ❖ When a molecule absorbs light having sufficient energy (e.g. UV-Vis radiation) to cause an electronic transitions, additional vibration and rotation transitions also occur
- ❖ Molecule can absorb one photon of just the right energy to cause the following simultaneous changes:



1. A transition from the ground electronic state E_0 to the E_1 excited electronic state
2. A change in the vibrational energy from the ground vibrational state of E_0 to an excited vibrational state of E_1
3. A transition from one rotational state of E_0 to a different rotational state of E_1
4. All the above transitions are quantized which means that they required certain exact amount of energy
5. Thus, total energy absorbed = $E_{\text{elec}} + E_{\text{vib}} + E_{\text{rot}}$
 $\Delta E_{\text{elec}} \gg \Delta E_{\text{vib}} \gg \Delta E_{\text{rot}}$


As a result, a large numbers of photons of certain wavelengths are absorbed by a molecule. These individual wavelengths are too numerous and too close to each other and a **spectrum of broad bands of absorbed wavelengths** are obtained




Spectrum (a graph that shows how absorbance varies with wavelength)

ELECTRONIC TRANSITIONS

POSSIBLE ELECTRONIC TRANSITIONS



IMPORTANT ELECTRONIC TRANSITIONS



CHROMOPHORE

The part of a molecule responsible for imparting color, are called as chromophores.

CHROMOPHORE

The functional groups containing multiple bonds capable of absorbing radiations.

e.g. NO_2 , $\text{N}=\text{O}$, $\text{C}=\text{O}$, $\text{C}=\text{N}$, $\text{C}\equiv\text{N}$, $\text{C}=\text{C}$, $\text{C}=\text{S}$

CHROMOPHORE

Chromophore	Example	Excitation	λ_{max} , nm	ϵ	Solvent
C=C	Ethene	$\pi \rightarrow \pi^*$	171	15,000	hexane
C \equiv C	1-Hexyne	$\pi \rightarrow \pi^*$	180	10,000	hexane
C=O	Ethanal	$n \rightarrow \pi^*$ $\pi \rightarrow \pi^*$	290 180	15 10,000	hexane hexane
N=O	Nitromethane	$n \rightarrow \pi^*$ $\pi \rightarrow \pi^*$	275 200	17 5,000	ethanol ethanol
C-X X=Br X=I	Methyl bromide Methyl iodide	$n \rightarrow \sigma^*$ $n \rightarrow \sigma^*$	205 255	200 360	hexane hexane

AUXOCHROME


The functional groups attached to a chromophore which modifies the ability of the chromophore to absorb light , altering the wavelength or intensity of absorption.

OR


The functional group with non-bonding electrons that does not absorb radiation in near UV region but when attached to a chromophore alters the wavelength & intensity of absorption.

AUXOCHROME


Benzene $\lambda_{\text{max}} = 255 \text{ nm}$





Phenol $\lambda_{\text{max}} = 270 \text{ nm}$



Aniline $\lambda_{\text{max}} = 280 \text{ nm}$









EFFECT OF CONJUGATION


Effect of Conjugation



Ethene




1,3-Butadiene





1,3,5-Hexatriene

Molecule	λ_{max} (nm)
$\text{CH}_2=\text{CH}_2$	171
$\text{CH}_2=\text{CH}-\text{CH}=\text{CH}_2$	217
$\text{CH}_2=\text{CH}-\text{CH}=\text{CH}-\text{CH}=\text{CH}_2$	258


Effect of Conjugation




Beta-Carotene



If a substance absorbs here...





B- Carotene



Lycopene

ABSORPTION & INTENSITY SHIFTS

- 1 • Bathochromic Shift (Red Shift)
- 2 • Hypsochromic Shift (Blue Shift)
- 3 • Hyperchromic Effect
- 4 • Hypochromic Effect

1

• Bathochromic Shift (Red Shift)

- When absorption maxima (λ_{max}) of a compound shifts to longer wavelength, it is known as bathochromic shift or red shift.
- The effect is due to presence of an auxochrome or by the change of solvent.
- e.g. An auxochrome group like $-\text{OH}$, $-\text{OCH}_3$ causes absorption of compound at longer wavelength.

2


• Hypsochromic Shift (Blue Shift)

- When absorption maxima (λ_{max}) of a compound shifts to shorter wavelength, it is known as hypsochromic shift or blue shift.
- The effect is due to presence of an group causes removal of conjugation or by the change of solvent.

3

• Hyperchromic Effect


- When absorption intensity (ϵ) of a compound is increased, it is known as hyperchromic shift.
- If auxochrome introduces to the compound, the intensity of absorption increases.



Pyridine

$$\lambda_{\max} = 257 \text{ nm}$$

$$\epsilon = 2750$$



2-methyl pyridine

$$\lambda_{\max} = 260 \text{ nm}$$


$$\epsilon = 3560$$

4

• Hypochromic Effect

- When absorption intensity (ϵ) of a compound is decreased, it is known as hypochromic shift.


Shifts and Effects



ULTRAVIOLET-VISIBLE
SPECTROSCOPY OF INORGANIC
COMPOUNDS

COMPOUNDS OF TRANSITION METAL COMPLEXES SOLUTION






Spectrochemical Series

$\text{I}^- < \text{Br}^- < \text{S}^{2-} < \text{SCN}^- < \text{Cl}^- < \text{N}_3^-$, $\text{F}^- < \text{urea}$, $\text{OH}^- < \text{ox}$, $\text{O}^{2-} < \text{H}_2\text{O} < \text{NCS}^- < \text{py}$, $\text{NH}_3 < \text{en} < \text{bpy}$, $\text{phen} < \text{NO}_2^- < \text{CH}_3^-$, $\text{C}_6\text{H}_5^- < \text{CN}^- < \text{CO}$.



An artist's wheel

Copyright © The McGraw-Hill Companies, Inc. Permission required for reproduction or display.





Absorbed color and observed colors

Absorbed Color	λ (nm)	Observed Color	λ (nm)
Violet	400	Green-yellow	560
Blue	450	Yellow	600
Blue-green	490	Red	620
Yellow-green	570	Violet	410
Yellow	580	Dark blue	430
Orange	600	Blue	450
Red	650	Green	520







FIVE BASIC OPTICAL INSTRUMENT COMPONENTS

- 1) **Source** - A stable source of radiant energy at the desired wavelength (or λ range).
- 2) **Sample Holder** - A transparent container used to hold the sample (cells, cuvettes, etc.).
- 3) **Wavelength Selector** - A device that isolates a restricted region of the EM spectrum used for measurement (monochromators, prisms, & filters).
- 4) **Detector** - Converts the radiant energy into a useable signal (usually electrical).
- 5) **Signal Processor & Readout** - Amplifies or attenuates the transduced signal and sends it to a readout device such as a meter, digital readout, chart recorder, computer, etc.

